20.7.2015

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Samples

1. Water (1L sample) from each heater-cooler unit. In particular, samples should include the filtered water supplies used to fill each/all tanks that are part of the unit, e.g. tanks for the cardioplegia fluid cooling circuit, tanks for the cardioplegia fluid heating circuits and tanks for the patient circuits.
2. Disinfectant solutions (50 ml) used for maintenance of the water tanks of the heater cooler-units.
3. Air (250-500 L samples) collected in the operating room during cardiovascular surgical interventions with CPB and during operation of the heater-cooling units. Two sampling positions at 0.3 m from the front and at 0.3 m behind the heater-cooler unit are recommended.

The air samples may be collected during a later phase of the investigation, targeting only those centres where *M. chimaera* has been detected in water samples from the heater-cooler units.

Culture methods

The above environmental water, aerosol and solution specimens should be processed in a biosafety containment level 2 mycobacteriology laboratory as described previously [11,23,24]. Whilst NTMs are not classified as Biosafety containment level 3 pathogens, given that mycobacteriology laboratories are of this standard (required for working with *M. tuberculosis*), further culturing and characterisation of samples may be processed in such a laboratory as this is where the laboratory resources and expertise is located.

Concentration of water samples can be performed by filtration or centrifugation as described previously [23,25]. The filtration method has provided a more sensitive detection of *M. chimaera* than centrifugation (Jakko van Ingen, personal communication).


In brief, environmental water samples can be processed as schematically described in the Annex-Figure1. If available, specimen decontamination of other environmental micro-organisms by use of cetylpyridinium chloride (0.005%) solution is to be preferred to the use of the more widely available standard solution of N-acetyl-L-cysteine sodium hydroxide (NALC-NaOH) or NaOH, as the latter decontamination solutions have been shown to decrease the recovery of MAC from environmental samples [25]. Water samples can be filtered and plated on selective (antimicrobial agent-containing) Middlebrook 7H11 or 7H10 agar or egg-based medium such as Löwenstein-Jensen agar. In addition, culture in liquid medium can be performed in parallel, using 50 ml samples centrifuged and then inoculated in the MGIT 960 system. Solid media should be preferred for culture if resources are limited.

Water samples should be collected preferably before the device is undergoing disinfection. Sterile sodium thiosulfate (20 mg/L) can be added to the water collection vials to neutralise any remaining hypochlorite in the water sample prior to further processing. However, this neutralisation step was not found to be necessary for recovery of MAC from potable water [24].

Processing of environmental air samples of bioaerosol can be performed as described previously [11,20,26,27].

Bioaerosol can be cultured from air samples collected with a surface impaction sampler. The sampler should be calibrated to sample a certain volume of air in a given period of time and utilise standard-sized selective (antimicrobial agent-containing) Middlebrook 7H11 or 7H10 agar plates. Samples of at least 250-500 L of air are recommended.

Figure 1. Methodology for environmental water sample processing and culture for MAC/ *Mycobacterium chimaera*

Abbreviations

CPC:	Cetylpyridinium chloride
SDW	Sterile distilled water
M7H10	Selective Middlebrook 7H10 agar
M7H11	Selective Middlebrook 7H11 agar
NTM	Non-tuberculous mycobacteria

References (excerpt)

11. Sax H, Bloemberg G, Hasse B, Sommerstein R, Kohler P, Achermann Y, et al. Prolonged outbreak of *Mycobacterium chimaera* infection after open chest heart surgery. Clin Infect Dis 2015; 61:67-75.
 20. Jagielski T, van Ingen J, Rastogi N, Dziadek J, Mazur PK, Bielecki J. Current methods in the molecular typing of *Mycobacterium tuberculosis* and other mycobacteria. Biomed Res Int 2014;645802.
 23. Parashar D, Chauhan DS, Sharma VD, Chauhan A, Chauhan SV, Katoch VM. Optimization of procedures for isolation of mycobacteria from soil and water samples obtained in northern India. Appl Environ Microbiol 2004; 70:3751-3.
 24. Thomson R, Carter R, Gilpin C, Coulter C, Hargreaves M. Comparison of methods for processing drinking water samples for the isolation of *Mycobacterium avium* and *Mycobacterium intracellulare*. Appl Environ Microbiol 2008; 74:3094-8.
 25. Radomski N, Cambau E, Moulin L, Haenn S, Moillon R, Lucas FS. Comparison of culture methods for isolation of nontuberculous mycobacteria from surface waters. Appl Environ Microbiol 2010; 76:3514-20.
 26. Public Health England. Protocol for environmental sampling, processing and culturing of water and air samples for the isolation of slow growing mycobacteria, June 2015. Available from: https://www.gov.uk/government/uploads/system/uploads/attachment_data/file/434197/Air_water_environmental_sampling_SoP.pdf.
 27. Zhang Y, Rajagopalan M, Brown BA, Wallace RJ, Jr. Randomly amplified polymorphic DNA PCR for comparison of *Mycobacterium abscessus* strains from nosocomial outbreaks. J Clin Microbiol 1997; 35:3132-9.
- New:** Kohler P, Kuster SP, Bloemberg G, Schulthess B, Frank M, Tanner FC, et al. Healthcare-associated prosthetic heart valve, aortic vascular graft, and disseminated *Mycobacterium chimaera* infections subsequent to open heart surgery. Eur Heart J 2015 Jul 17. doi:10.1093/eurheartj/ehv342
-